

HAGGAI
SURGICAL

ELECTROSURGICAL

INSTRUMENTS

CATALOGUE



Table of Contents

❑ European Bipolar Forceps.....	1-6
❑ European Irrigation Bipolar Forceps.....	7-7
❑ European Non-Stick Bipolar Forceps.....	8-10
❑ US Bipolar Forceps.....	11-16
❑ Monopolar Forceps.....	17-17
❑ US Irrigation Bipolar Forceps.....	18-18
❑ US Non-Stick Bipolar Forceps.....	19-20
❑ Hand-Switch Bipolar Forceps.....	21-21
❑ Hand-Switch Monopolar Forceps.....	22-22
❑ Diathermy Instruments.....	23-24
❑ Disposable Bipolar Forceps.....	25-26
❑ OB/Gyn Leep Instruments.....	27-29
❑ OB/Gyn Laser Surgery Instruments.....	30-31
❑ Laparoscopic & Arthroscopic Electrode.....	32-33
❑ Electrosurgical Pencil & Electrode.....	34-38
❑ Electrosurgical Cable.....	39-41
❑ Bipolar Artery Sealer.....	42-42



The haggai surgical forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.

Advantages of Reusable European Bipolar Forceps

- ✓ Highly polished tips reduce tissue adherence.
- ✓ Light closing force causes less fatigue over time.
- ✓ Nylon coating resists abrasion and provides superb insulation.
- ✓ Precision tip alignment aids fine control.
- ✓ Reusable, durable, and designed for multiple uses.
- ✓ The Instruments is autoclavable.
- ✓ Compatible with all common electrosurgery unit.

02 EUROPEAN BIPOLAR FORCEPS



HS-381-101 McPherson



HS-381-102 McPherson Ang.



HS-381-104 Jeweller #5 Str



HS-381-106 Jeweller #5 Ang




HS-381-106 Jeweller #5 Ang




HS-381-110 Adson

Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-100	McPherson	0.5	10.2
HS-381-102	McPherson Angled	0.5	10.2
HS-381-104	Jeweller #5	0.5, 0.7	11.5
HS-381-106	Jeweller #5 Angled	0.5, 0.7	11.5
HS-381-108	Jeweller #7	0.5, 0.7	11.5
HS-381-110	Adson	1.0	12.0
HS-381-112	Adson	1.0	15.0



Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref		Tip Size (mm)	Length Size (cm)
HS-381-114	Straight	0.5, 1.0, 2.0	15.0
HS-381-116	Straight	0.5, 1.0, 2.0	18.0
HS-381-118	Straight	0.5, 1.0, 2.0	20.0
HS-381-120	Straight	0.5, 1.0, 2.0	22.0
HS-381-122	Straight	0.5, 1.0, 2.0	24.0




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-124	Straight Angled	0.5, 1.0, 2.0	15.0
HS-381-126	Straight Angled	0.5, 1.0, 2.0	18.0
HS-381-128	Straight Angled	0.5, 1.0, 2.0	20.0
HS-381-130	Straight Angled	0.5, 1.0, 2.0	22.0
HS-381-132	Straight Angled	0.5, 1.0, 2.0	24.0




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-134	Bayonet	0.5, 1.0, 2.0	16.5
HS-381-136	Bayonet	0.5, 1.0, 2.0	18.0
HS-381-138	Bayonet	0.5, 1.0, 2.0	20.0
HS-381-140	Bayonet	0.5, 1.0, 2.0	22.0
HS-381-142	Bayonet	0.5, 1.0, 2.0	24.0



Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-144	Bayonet Angled Up	0.5, 1.0, 2.0	16.5
HS-381-146	Bayonet Angled Up	0.5, 1.0, 2.0	18.0
HS-381-148	Bayonet Angled Up	0.5, 1.0, 2.0	20.0
HS-381-150	Bayonet Angled Up	0.5, 1.0, 2.0	22.0
HS-381-152	Bayonet Angled Up	0.5, 1.0, 2.0	24.0



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size
HS-381-154	Bayonet Angled Down	0.5, 1.0, 2.0	16.5
HS-381-156	Bayonet Angled Down	0.5, 1.0, 2.0	18.0
HS-381-158	Bayonet Angled Down	0.5, 1.0, 2.0	20.0
HS-381-160	Bayonet Angled Down	0.5, 1.0, 2.0	22.0
HS-381-162	Bayonet Angled Down	0.5, 1.0, 2.0	24.0



It is well known that depending on the generator and the surface quality of the forceps jaws, coagulated tissue tends to more or less stick to the instrument jaw. This unwanted effect can occur particularly often with liver tissue and is critical in neurosurgery. For these applications in particular, we have designed our bipolar forceps with integrated irrigation channel. The coagulation area is irrigated through a built-in irrigation channel inside the forceps branch which can be connected to a standard irrigation hose. By means of this irrigation the sticking of the tissue is substantially reduced and surgical work is made easier.

Advantages of Reusable European Irrigation Bipolar Forceps

- ✓ A universal hub connector allows irrigation tubing to be connected easily.
- ✓ The ergonomically designed handle provides less fatigue and a well-balanced instrument.
- ✓ Nylon coating resists abrasion and provides superb insulation.
- ✓ For high visibility, our irrigation tube can be recessed into the forceps blank (buried irrigation).
- ✓ Compatible with all common electrosurgery units.


HS-381-170-IR Bayonet Irrigation

HS-381-168-IR Bayonet Irrigation

HS-381-166-IR Straight Irrigation

HS-381-164-IR Straight Irrigation

Tip Chart



Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-164-IR	Straight Irrigation	0.5, 1.0, 2.0	18.0
HS-381-166-IR	Straight Irrigation	0.5, 1.0, 2.0	20.0
HS-381-168-IR	Bayonet Irrigation	0.5, 1.0, 2.0	18.0
HS-381-170-IR	Bayonet Irrigation	0.5, 1.0, 2.0	20.0

Note: More Bipolar Irrigation Forceps are available use IR extension to any regular bipolar forceps.

Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm




During the caogulation process, tissue and tips of forceps can stick together. Haggai non-stick bipolar forceps are made with a special silver alloy and therefore prevent tissue from sticking. Disturbing and time consuming disruptions for cleaning the tips during surgery are thus avoided. Due to the good thermal qualities of silver the temperature of the forceps tips does not exceed 100°C during coagulation, which enhances the non-stick effect. With over ten years of experience we can guarantee that our bipolar forceps do not stick. The non-stick effect remains constant and is not reduced by frequent sterilization. Our non-stick bipolar forceps together with the appropriate connection cables are compatible with all common types of high-frequency units such as Erbe, Martin, Berchtold, Valleylab. Crosstex Bipolar forceps non-stick are available in straight, angled or bayonet shape. Available lengths range from 110 to 300 mm. Tips can be made from 0.25 to 2.0 mm width. Forceps are also available with different measurements as well as with flat pin or 2-pin plug.


Advantages of Reusable Non-Stick European Bipolar Forceps

- ✓ Ten years of experience speak for themselves. Our forceps do not stick – never!
- ✓ No cleaning during surgery
- ✓ Especially high thermal qualities
- ✓ All Non-stick Bipolar Forceps are re-usable and autoclavable.
- ✓ Compatible with all common electrosurgery unit.

11 EUROPEAN NON-STICK BIPOLAR FORCEP




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-172-NS	McPherson	0.5,	10.2
HS-381-174-NS	McPherson Angled	0.5,	10.2
HS-381-176-NS	Jeweller #5	0.5,0.7	11.5
HS-381-178-NS	Jeweller #5 Angled	0.5,0.7	11.5
HS-381-180-NS	Adson	1.0	12.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-182-NS	Straight	0.5, 1.0, 2.0	15.0
HS-381-184-NS	Straight	0.5, 1.0, 2.0	18.0
HS-381-186-NS	Straight	0.5, 1.0, 2.0	20.0
HS-381-188-NS	Straight	0.5, 1.0, 2.0	22.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-190-NS	Bayonet	0.5, 1.0, 2.0	16.5
HS-381-192-NS	Bayonet	0.5, 1.0, 2.0	18.0
HS-381-194-NS	Bayonet	0.5, 1.0, 2.0	20.0
HS-381-196-NS	Bayonet	0.5, 1.0, 2.0	22.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.




The HAGGAI SURGICAL forceps are available in various lengths and shapes for different surgical and anatomical requirements. The reusable Bipolar forceps make it easy to coagulate tissue. Sticking effects on the tips of the forceps are reduced to a minimum.


Advantages of Reusable US OPHTHALMIC Bipolar Forceps

- ✓ Highly polished tips reduce tissue adherence.
- ✓ Light closing force causes less fatigue over time.
- ✓ Nylon coating resists abrasion and provides superb insulation.
- ✓ Reusable, durable, and designed for multiple uses.
- ✓ Precision tip alignment aids fine control.
- ✓ The Instruments is autoclavable.
- ✓ Compatible with all common electrosurgery unit.

15 US OPHTHALMIC BIPOLAR FORCEPS




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-230	McPherson	0.5	10.2
HS-381-232	McPherson Angled	0.5	10.2
HS-381-234	Coaptation Nedler	0.5	10.2
HS-381-236	Coaptation	0.5	10.8
HS-381-238	Jeweller #5	Needle, 0.5	11.5

16 US OPHTHALMIC BIPOLAR FORCEPS



Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-240	Jeweller #4	Needle, 0.5	11.5
HS-381-242	Tenzel	Needle, 0.5	10.2
HS-381-244	Jeweller #7	Needle, 0.5	11.5
HS-381-246	Iris Straight	Needle, 0.5	10.8
HS-381-248	Iris Angled	0.5, 1.0	10.8




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-250	Adson	1.0	12.0
HS-381-252	Semkin	0.5, 1.0	14.0
HS-381-254	Cushing	1.0, 1.5, 2.0	15.9
HS-381-256	Cushing	1.0, 1.5, 2.0	17.8




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-258	Cushing Gursh	1.0, 1.5, 2.0	17.8
HS-381-260	Gerald Straight	1.0, 1.5, 2.0	17.8
HS-381-262	Gerald Angled	1.0, 1.5, 2.0	17.8
HS-381-264	Gerald Bayonet	1.0, 1.5, 2.0	18.4




Tip Chart




Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-266	Jansen Bayonet	1.0, 1.5, 2.0	15.2
HS-381-268	Jansen Bayonet	1.0, 1.5, 2.0	17.8
HS-381-270	Jansen Bayonet	1.0, 1.5, 2.0	20.4
HS-381-272	Malis	1.0, 1.5, 2.0	17.8



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #:
TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-274	Scoville-Greenwood	1.0, 1.5	19.1
HS-381-276	Yasargil	1.0, 1.5	19.7
HS-381-278	Hardy Bayonet	1.0, 1.5	21.6
HS-381-280	Hardy Bayonet Angled	1.0, 1.5	21.6



HAGGAI SURGICAL Monopolar instruments are tested at high volts to ensure proper functioning during the most stringent operations. With the best quality stainless steel, excellent finishing and superior insulation, we provide worldwide guarantee for all our Monopolar instruments.

Advantages of Reusable Monopolar Forceps

- ✓ Lightweight and well balanced to prevent fatigue.
- ✓ Highly polished tips minimise sticking.
- ✓ Nylon coating resists abrasion and provides superb insulation.
- ✓ Precision tip alignment aids fine control.
- ✓ Compatible with all common electrosurgery unit.


HS-381-282 Jeweller #4



HS-381-284 Jeweller #7



HS-381-286 Adson



HS-381-288 Potts-Smith



HS-381-292 Potts-Smith



HS-381-294 Jansen Bayonet




HS-381-296 Jansen Bayonet



HS-381-302 Jansen Bayonet UP



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-282	Jeweller #4	Needle	11.5
HS-381-284	Jeweller #7	Needle	11.5
HS-381-286	Adson	1.0	12.0
HS-381-288	Potts-Smith	2.0	17.8
HS-381-290	Potts-Smith	2.0	21.0
HS-381-292	Potts-Smith	2.0	24.1
HS-381-294	Jansen Bayonet	1.5	13.3
HS-381-296	Jansen Bayonet	1.5	15.9
HS-381-298	Jansen Bayonet	2.0	19.1
HS-381-300	Jansen Bayonet	2.0	21.0
HS-381-302	Jansen Bayonet UP	2.0	19.1
HS-381-304	Jansen Bayonet UP	2.0	21.0



It is well known that depending on the generator and the surface quality of the forceps jaws, coagulated tissue tends to more or less stick to the instrument jaw. This unwanted effect can occur particularly often with liver tissue and is critical in neurosurgery. For these applications in particular, we have designed our Bipolar Forceps with integrated irrigation channel. The coagulation area is irrigated through a built-in irrigation channel inside the forceps branch which can be connected to a standard irrigation hose. By means of this irrigation the sticking of the tissue is substantially reduced and surgical work is made easier.

Advantages of Reusable US Irrigation Forceps

- ✓ A universal hub connector allows irrigation tubing to be connected easily.
- ✓ The ergonomically designed handle provides less fatigue and a well balanced instrument.
- ✓ For high visibility our irrigation tube can be recessed into the forceps blank (buried irrigation).
- ✓ Compatible with all common electrosurgery unit.



HS-381-306-IR Hardy Bayonet



HS-381-308-IR Jansen Bayonet



HS-381-310-IR Scoville-Greenwood




HS-381-312-IR Gerald Bayonet



HS-381-314-IR Straight


Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-306-IR	Hardy Bayonet	0.5, 1.0, 1.5	21.6
HS-381-308-IR	Jansen Bayonet	0.5, 1.0, 1.5	19.1
HS-381-310-IR	Scoville-Greenwood	0.5, 1.0, 1.5	19.7
HS-381-312-IR	Gerald Bayonet	0.5, 1.0, 1.5	18.4
HS-381-314-IR	Straight	0.5, 1.0, 1.5	15.9

Note: More Irrigation Bipolar Forceps are available use IR extension to any regular bipolar forceps.



During the coagulation process, tissue and tips of forceps can stick together. Haggai Surgical non-stick bipolar forceps are made with a special silver alloy and therefore prevent tissue from sticking. Disturbing and time consuming disruptions for cleaning the tips during surgery are thus avoided. Due to the good thermal qualities of silver the temperature of the forceps tips does not exceed 100°C during coagulation, which enhances the non-stick effect. With over ten years of experience we can guarantee that our bipolar forceps do not stick. The non-stick effect remains constant and is not reduced by frequent sterilization. Our non-stick bipolar forceps together with the appropriate connection cables are compatible with all common types of high-frequency units such as Erbe, Martin, Haggai surgical told, Valleylab. Haggai Surgical Bipolar forceps non-stick are available in straight, angled or bayonet shape. Available lengths range from 110 to 300 mm. Tips can be made from 0.25 to 2.0 mm width. Forceps are also available with different measurements as well as with flat pin or 2-pin plug.

Advantages of Reusable Non-Stick European Bipolar Forceps

- ✓ Ten years of experience speak for themselves. Our forceps do not stick – never!
- ✓ No cleaning during surgery.
- ✓ Especially high thermal qualities.
- ✓ All Non-stick Bipolar Forceps are re-usable and autoclavable.
- ✓ *Compatible with all common electrosurgery unit.*



HS-381-316-NS McPherson



HS-381-318-NS McPherson Angled



HS-381-320-NS Jeweller #5



HS-381-322-NS Jeweller #4




HS-381-324-NS Adson



HS-381-326-NS Semkin

Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-316-NS	McPherson	0.5	10.2
HS-381-318-NS	McPherson Angled	0.5	10.2
HS-381-320-NS	Jeweller #5	Needle, 0.5	11.5
HS-381-322-NS	Jeweller #4	Needle, 0.5	11.5
HS-381-324-NS	Adson	1.0	12.0
HS-381-326-NS	Semkin	0.5, 1.0	14.0

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.



HS-381-328-NS Cushing



HS-381-332-NS Gerald Bayonet



HS-381-334-NS Jansen Bayonet



HS-381-336-NS Scoville-Greenwood




HS-381-338-NS Yasargil



HS-381-340-NS Hardy Bayonet

Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-328-NS	Cushing	1.0, 1.5, 2.0	15.9
HS-381-330-NS	Cushing	1.0, 1.5, 2.0	17.8
HS-381-332-NS	Gerald Bayonet	1.0, 1.5	18.4
HS-381-334-NS	Jansen Bayonet	1.0, 1.5	15.9
HS-381-336-NS	Scoville-Greenwood	1.0, 1.5	19.1
HS-381-338-NS	Yasargil	1.0, 1.5	19.7
HS-381-340-NS	Hardy Bayonet	1.0, 1.5	21.6

Note: More Non-Stick Bipolar Forceps are available use NS extension to any regular bipolar forceps.

HS-381-341-HS Jeweller #4


HS-381-342-HS Semkin

HS-381-343-HS Jansen Bayonet

HS-381-344-HS Scoville-Greenwood



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm


Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-341-HS	Jeweller #4	0.5,	11.5
HS-381-342-HS	Semkin	0.7	14.0
HS-381-343-HS	Jansen Bayonet	0.7	19.1
HS-381-344-HS	Scoville-Greenwood	1.5	19.7

HS-381-345-HS Jansen Bayonet


HS-381-346-HS Jansen Bayonet

HS-381-347-HS Potts-Smith

HS-381-348-HS Potts-Smith



Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-345-HS	Jansen Bayonet	2.0	19.0
HS-381-346-HS	Jansen Bayonet	2.0	21.0
HS-381-347-HS	Potts-Smith	2.0	21.0
HS-381-348-HS	Potts-Smith	2.0	24.0



Our Diathermy instruments are tested at high volts to ensure proper functioning during the most stringent operations. With the best quality stainless steel, excellent finishing and superior insulation, we provide worldwide guarantee for all our Diathermy instruments. Diathermy instruments are available in four base types. Models popular Diathermy Instruments : Adson, Jeweller, McIndoe, Standard, Bayonet, Artery Forceps, Wilson Hey Forceps, Diathermy Scissors, etc in various lengths, tip sizes, and 4 plug types

Advantages of Reusable Diathermy Instruments

- ✓ Lightweight and well balanced to prevent fatigue.
- ✓ Highly polished tips minimise sticking.
- ✓ Precision tip alignment aids fine control.
- ✓ Compatible with all common electro surgery unit.



Tuner Warwick Scissors (STR)
HS-381-350



Tuner Warwick Scissors (CRV)
HS-381-352



Dressing Scissors (STR)
HS-381-358



Dressing Scissors (CRV)
HS-381-360



Wilson Hay Forceps (STR)
HS-381-366



Wilson Hay Forceps (CRV)
HS-381-368



Artery Forceps (STR)
HS-381-374



Artery Forceps (CRV)
HS-381-376



Abbey Needle
HS-381-382



Richie Forceps
HS-381-384

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-350	Tuner Warwick Scissors (STR)	-	17.8
HS-381-352	Tuner Warwick Scissors (CRV)	-	17.8
HS-381-354	Tuner Warwick Scissors (STR)	-	23.0
HS-381-356	Tuner Warwick Scissors (CRV)	-	23.0
HS-381-358	Dressing Scissors (STR)	-	17.8
HS-381-360	Dressing Scissors (CRV)	-	17.8
HS-381-362	Dressing Scissors (STR)	-	23.0
HS-381-364	Dressing Scissors (CRV)	-	23.0
HS-381-366	Wilson Hay Forceps (STR)	-	17.8
HS-381-368	Wilson Hay Forceps (CRV)	-	17.8
HS-381-370	Wilson Hay Forceps (STR)	-	23.0
HS-381-372	Wilson Hay Forceps (CRV)	-	23.0
HS-381-374	Artery Forceps (STR)	-	17.8
HS-381-376	Artery Forceps (CRV)	-	17.8
HS-381-378	Artery Forceps (STR)	-	23.0
HS-381-380	Artery Forceps (CRV)	-	23.0
HS-381-382	Abbey Needle	-	17.5
HS-381-384	Richie Forceps	-	17.8
HS-381-386	Richie Forceps	-	20.0
HS-381-388	Richie Forceps	-	23.0

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-390	Tuner Warwick Scissors (STR)	-	17.8
HS-381-392	Tuner Warwick Scissors (CRV)	-	17.8
HS-381-394	Tuner Warwick Scissors (STR)	-	23.0
HS-381-396	Tuner Warwick Scissors (CRV)	-	23.0
HS-381-398	Dressing Scissors (STR)	-	17.8
HS-381-400	Dressing Scissors (CRV)	-	17.8
HS-381-402	Dressing Scissors (STR)	-	23.0
HS-381-404	Dressing Scissors (CRV)	-	23.0
HS-381-406	Wilson Hay Forceps (STR)	-	17.8
HS-381-408	Wilson Hay Forceps (CRV)	-	17.8
HS-381-410	Wilson Hay Forceps (STR)	-	23.0
HS-381-412	Wilson Hay Forceps (CRV)	-	23.0
HS-381-414	Artery Forceps (STR)	-	17.8
HS-381-416	Artery Forceps (CRV)	-	17.8
HS-381-418	Artery Forceps (STR)	-	23.0
HS-381-420	Artery Forceps (CRV)	-	23.0
HS-381-422	Abbey Needle	-	17.5
HS-381-424	Richie Forceps	-	17.8
HS-381-426	Richie Forceps	-	20.0
HS-381-428	Richie Forceps	-	23.0



Adson
HS-381-430



McINDOE
HS-381-432



Gillies
HS-381-434



Dissecting
HS-381-442



Gerald
HS-381-450



Adson Bayonet
HS-381-452



Adson Bayonet UP
HS-381-456



Remington Hobb
HS-381-462

ENGLISH MONOPOLAR MALE 4.8MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
HS-381-430	Adson	1.0	12.0
HS-381-432	McINDOE	1.0,2.0	15.0
HS-381-434	Gillies	-	15.0
HS-381-436	Gillies	-	17.8
HS-381-438	Gillies	-	20.0
HS-381-440	Gillies	-	24.0
HS-381-442	Dissecting	1.0,2.0	15.0
HS-381-444	Dissecting	1.0,2.0	17.8
HS-381-446	Dissecting	1.0,2.0	21.0
HS-381-448	Dissecting	1.0,2.0	24.1
HS-381-450	Gerald 1x2	-	17.8
HS-381-452	Adson Bayonet	1.0,2.0	16.0
HS-381-454	Adson Bayonet	1.0,2.0	19.1
HS-381-456	Adson Bayonet Angled UP	1.0,2.0	16.0
HS-381-458	Adson Bayonet Angled UP	1.0,2.0	19.1
HS-381-460	Adson Bayonet Angled UP	1.0,2.0	21.0
HS-381-462	Remington Hobb	1.0,2.0	24.1

EUROPEAN MONOPOLAR FEMALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
HS-381-364	Adson	1.0	12.0
HS-381-366	McINDOE	1.0,2.0	15.0
HS-381-368	Gillies	-	15.0
HS-381-370	Gillies	-	17.8
HS-381-372	Gillies	-	20.0
HS-381-474	Gillies	-	24.0
HS-381-476	Dissecting	1.0,2.0	15.0
HS-381-478	Dissecting	1.0,2.0	17.8
HS-381-480	Dissecting	1.0,2.0	21.0
HS-381-482	Dissecting	1.0,2.0	24.1
HS-381-484	Gerald 1x2	-	17.8
HS-381-486	Adson Bayonet	1.0,2.0	16.0
HS-381-488	Adson Bayonet	1.0,2.0	19.1
HS-381-490	Adson Bayonet Angled UP	1.0,2.0	16.0
HS-381-492	Adson Bayonet Angled UP	1.0,2.0	19.1
HS-381-494	Adson Bayonet Angled UP	1.0,2.0	21.0
HS-381-496	Remington Hobb	1.0,2.0	24.1

ENGLISH MONOPOLAR MALE 4.0MM CONNECTOR

Ref	Item Name	Tip Size	Length Size
HS-381-498	Adson	1.0	12.0
HS-381-500	McINDOE	1.0,2.0	15.0
HS-381-502	Gillies	-	15.0
HS-381-504	Gillies	-	17.8
HS-381-506	Gillies	-	20.0
HS-381-508	Gillies	-	24.0
HS-381-510	Dissecting	1.0,2.0	15.0
HS-381-512	Dissecting	1.0,2.0	17.8
HS-381-514	Dissecting	1.0,2.0	21.0
HS-381-516	Dissecting	1.0,2.0	24.1
HS-381-518	Gerald 1x2	-	17.8
HS-381-520	Adson Bayonet	1.0,2.0	16.0
HS-381-522	Adson Bayonet	1.0,2.0	19.1
HS-381-524	Adson Bayonet Angled UP	1.0,2.0	16.0
HS-381-526	Adson Bayonet Angled UP	1.0,2.0	19.1
HS-381-528	Adson Bayonet Angled UP	1.0,2.0	21.0
HS-381-530	Remington Hobb	1.0,2.0	24.1

MONOPOLAR WITHOUT CONNECTOR

Ref	Item Name	Tip Size	Length Size
HS-381-532	Adson	1.0	12.0
HS-381-534	McINDOE	1.0,2.0	15.0
HS-381-536	Gillies	-	15.0
HS-381-5HS-38	Gillies	-	17.8
HS-381-540	Gillies	-	20.0
HS-381-542	Gillies	-	24.0
HS-381-544	Dissecting	1.0,2.0	15.0
HS-381-546	Dissecting	1.0,2.0	17.8
HS-381-548	Dissecting	1.0,2.0	21.0
HS-381-550	Dissecting	1.0,2.0	24.1
HS-381-552	Gerald 1x2	-	17.8
HS-381-554	Adson Bayonet	1.0,2.0	16.0
HS-381-556	Adson Bayonet	1.0,2.0	19.1
HS-381-558	Adson Bayonet Angled UP	1.0,2.0	16.0
HS-381-560	Adson Bayonet Angled UP	1.0,2.0	19.1
HS-381-562	Adson Bayonet Angled UP	1.0,2.0	21.0
HS-381-564	Remington Hobb	1.0,2.0	24.1



Haggai surgical disposable bipolar forceps have been designed and tested to meet the practical. The main purpose for offering single-use bipolar forceps is to avoid the risk of contamination while remaining cost effective and utilizing the functionalities of the higher end forceps. Hardened stainless steel provides an optimum spring tension and serves as the basis for a superior instrument.

Advantages of Disposable Bipolar Forceps

- ✓ Lightweight and well balanced to prevent fatigue.
- ✓ Highly polished tips minimize sticking.
- ✓ Nylon coating resists abrasion and provides superb insulation.
- ✓ Precision tip alignment aids fine control.
- ✓ Cost effective in comparison to re-usable.
- ✓ All Bipolar forceps available with cable, sterile packed.
- ✓ Compatible with all common electro surgery unit.



HS-381-570-DP Jeweller #5



HS-381-572-DP Straight



HS-381-574-DP Bayonet




HS-381-576-DP Straight



HS-381-578-DP Straight Angled

Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-570-DP	Jeweller #5 With 3m Cable	0.5	11.5
HS-381-572-DP	Straight With 3m Cable	1.0, 1.5	18.0
HS-381-574-DP	Bayonet With 3m Cable	1.0, 1.5	18.0
HS-381-576-DP	Straight	1.0, 1.5	18.0
HS-381-578-DP	Bayonet	1.0, 1.5	20.0

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.



HS-381-580-DP
Gerald Straight


HS-381-582-DP
Jansen Irrigation

HS-381-584-DP
Jansen Bayonet

HS-381-586-DP
Jeweller #5

HS-381-588-DP
Adson

Tip Chart



Tip size can be changed as per your requirement, just add the following extensions fit tip size at the end of Ref #: TIP- Needle, 05, 0.7mm, 1.0mm, 1.5mm & 2.0mm

Ref	Item Name	Tip Size (mm)	Length Size (cm)
HS-381-580-DP	US. Gerald	1.0, 1.5	17.8
HS-381-582-DP	US. Jansen Bayonet Irrigation	1.0, 1.5	19.1
HS-381-584-DP	US. Jansen Bayonet	1.0, 1.5	17.8
HS-381-586-DP	US. Jeweller #5	0.5	11.5
HS-381-588-DP	US. Adson	1.0	12.1

Note: More Disposable Bipolar Forceps are Available Use DP extension to any regular Bipolar Forceps.



The first choice for proven performance, Haggai surgical provides the broadest selection of fully autoclavable instruments for LEEP procedures. Our "rose coated" instruments will provide you with durability and performance unmatched in the industry. Models Available: Graves, Pederson, VU-More, 4-Way, Cusco, Collin, Weisman, Gelpi, Rigby, Lateral, Lateral Vaginal, Endocervical, Iris, Emmet, Straight, two-prong, three-prong, etc.



Disposable Tube



Graves Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-600**

Medium **HS-381-602**

Large **HS-381-604**

Small

HS-381-606

Medium

HS-381-608

Large

HS-381-610



Disposable Tube



Pederson Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-612**

Medium **HS-381-614**

Large **HS-381-616**

Small

HS-381-618

Medium

HS-381-620

Large

HS-381-622



Disposable Tube



VU-Max Grave Vaginal Speculum

With Smoke Evacuation Tube

Medium **HS-381-624**

Large **HS-381-626**

Medium

HS-381-628

Large

HS-381-630



Disposable Tube



VU-Max Pederson Vaginal Speculum

With Smoke Evacuation Tube

Medium **HS-381-632**

Large **HS-381-634**

Medium

HS-381-636

Large

HS-381-638



Collin Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-640**

Medium **HS-381-642**

Large **HS-381-644**



CUSCO Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-646**

Medium **HS-381-648**

Large **HS-381-650**



VU-More Graves Vaginal Speculum
With Smoke Evacuation Tube
 Medium **HS-381-652**
 Large **HS-381-654**

Disposable Tube



Medium
HS-381-656
 Large
HS-381-658



4-Blades Speculum
With Smoke Evacuation Tube
 Medium **HS-381-660**
 Large **HS-381-662**



**Weisman Vaginal Speculum
 Right Side Open**
With Smoke Evacuation Tube
 Medium **HS-381-664**
 Large **HS-381-666**

Disposable Tube



Medium
HS-381-668
 Large
HS-381-670



**Weisman Vaginal Speculum
 Left Side Open**
With Smoke Evacuation Tube
 Medium **HS-381-672**
 Large **HS-381-674**

Disposable Tube



Medium
HS-381-676
 Large
HS-381-678



**Lateral Retractor
 Wide Open**
HS-381-680



**Rigby Retractor
 Wide Open**
HS-381-682



**Gelpi Retractor
 Wide Open**
HS-381-684



Kogan Speculum
 3mm **HS-381-686**
 4mm **HS-381-688**
 5mm **HS-381-690**



**Lateral Retractor
 21cm 22mmx75mm**
HS-381-692



**Schroeder
Tanaculum 24cm
HS-381-694**



**Schroeder
Tanaculum 24cm
HS-381-696**



**Single Tooth
Tanaculum 24cm
HS-381-698**



**Campion
Forceps 22cm
HS-381-700**



**Spong
Forceps 22cm
HS-381-702**

**Dressing Forceps HS-
381-704 21cm
HS-381-706 25cm**



**Tissue Forceps HS-
381-708 21cm
HS-381-710 25cm**



**Iris Hook
HS-381-712 25cm**



**Emmet Hook
HS-381-714 25cm**



**Emmet Hook Straight
HS-381-716 25cm**



**Emmet Hook Two Prong
HS-381-718 25cm**



**Emmet Hook Three Prong
HS-381-720 25cm**





Graves Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-750**

Medium **HS-381-752**

Large **HS-381-754**



Lateral Vaginal Speculum

HS-381-766



Pederson Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-756**

Medium **HS-381-758**

Large **HS-381-760**




Collin Vaginal Speculum

With Smoke Evacuation Tube

Small **HS-381-762**

Medium **HS-381-764**

Large **HS-381-766**



Kogan Endocervical Speculum



3mm
HS-381-770



4mm
HS-381-772



5mm
HS-381-774



Rigby Retractor
HS-381-778



Gelpi Retractor
HS-381-776



Schroeder Tenaculum
HS-381-780



Iris Hook
HS-381-782



Emmet Hook
HS-381-784




Emmet Hook Straight
HS-381-786



Emmet Hook Two Prong
HS-381-788



Emmet Hook Three Prong
HS-381-790



These E-Z Clean Laparoscopic Electrodes are made from very high quality raw material which ensures hassle free work performance at its user end. These E-Z Clean Laparoscopic Electrodes are widely finds its applications in various medical sectors.

Salient Features of our E-Z Clean Laparoscopic Electrodes

- ✓ *Provide excellent dissection and coagulation for a variety of surgical applications.*
- ✓ *High-tech PTFE coating for reduced eschar buildup.*
- ✓ *Reintroduction during surgery is minimized.*
- ✓ *Lower power settings for safer delivery of electrosurgery.*
- ✓ *Tips are cost-effective because they Save Time and Money.*

Without Suction



L-Hook Electrode 36cm
HS-381-800




J-Hook Electrode 36cm
HS-381-802



Needle Electrode 36cm
HS-381-804



Spectula Electrode 36cm
HS-381-806



Ball Electrode 36cm
HS-381-808



With Suction



L-Hook Electrode 36cm
HS-381-810



J-Hook Electrode 36cm
HS-381-812



Needle Electrode 36cm
HS-381-814



Spectula Electrode 36cm
HS-381-816



Ball Electrode 36cm
HS-381-818



Dagger 12xm
HS-381-820



Dagger 18cm
HS-381-822



Hook 12cm
HS-381-824



Hook 18cm
HS-381-826



Hockey 12cm
HS-381-828



Hockey 18cm
HS-381-830



Frazier Suction Tube 11FR 13.5cm
HS-381-832



Frazier Suction Tube 11FR 20.5cm
HS-381-834



Frazier Suction Tube 11FR 40.6cm
HS-381-836





Electrode handle with blade electrode for single-use. Sterile packaging in individual pouches. Ready to use. This handle is made of high quality material and fits into the surgeon's hand comfortably due to its ergonomic shape. The quality meets highest expectations and fulfills all required international safety regulations.

- ✓ Single-use electrode handle including blade electrode and 3 m cable (also available: 5 m)
- ✓ Individually packaged, sterile.
- ✓ Push buttons for cutting and coagulation.
- ✓ Unit side: 3-pin connector.
 - ✓ Tips are cost-effective because they Save Time and Money.
- ✓ Suitable also for use with ERBE, Martin or Berchtold units, by addition of an adapter especially designed for such purpose.



Rocker Switch Pencil (Single Use)
3m Cable HS-381-854



Foot Control Pencil (Single Use)
3m Cable HS-381-856
5m Cable HS-381-858



Finger Switch Pencil (Reusable)
3m Silicon Cable HS-381-860



Finger Switch Pencil (Single Use)
3m Cable HS-381-850
5m Cable HS-381-852



Electrode Handle (Reusable)
3m Silicon Cable HS-381-862



Blade Electrode
 Disposable: 381-900-S Reusable: 381-900-R



Blade Electrode Crv
 Disposable: 381-902-S Reusable: 381-902-R



Needle Electrode
 Disposable: 381-904-S Reusable: 381-904-R



Needle Electrode Crv
 Disposable: 381-906-S Reusable: 381-906-R



Ball Electrode 2mm
 Disposable: 381-908-S Reusable: 381-908-R



Ball Electrode Crv 2mm
 Disposable: 381-910-S Reusable: 381-910-R



Ball Electrode 3mm
 Disposable: 381-912-S Reusable: 381-912-R



Ball Electrode 3mm Crv
 Disposable: 381-914-S Reusable: 381-914-R



Ball Electrode 4mm
 Disposable: 381-916-S Reusable: 381-916-R



Ball Electrode 4mm Crv
 Disposable: 381-918-S Reusable: 381-918-R



Ball Electrode 5mm
 Disposable: 381-920-S Reusable: 381-920-R



Ball Electrode 5mm Crv
 Disposable: 381-922-S Reusable: 381-922-R



Blade Electrode Shaft 13cm
 Disposable: 381-924-S Reusable: 381-924-R



Needle Electrode Shaft 13cm
 Disposable: 381-926-S Reusable: 381-926-R



Ball Electrode 3mm Shaft 13cm
 Disposable: 381-928-S Reusable: 381-928-R



Ball Electrode 5mm Shaft 13cm
 Disposable: 381-930-S Reusable: 381-930-R



Conization Electrode
 Disposable: 381-932-S Reusable: 381-932-R



Circle Loop Electrode 4mm
 Disposable: 381-934-S Reusable: 381-934-R



Circle Loop Electrode 6mm
 Disposable: 381-936-S Reusable: 381-936-R



Circle Loop Electrode 8mm
 Disposable: 381-938-S Reusable: 381-938-R



Diamond Electrode
 Disposable: 381-940-S Reusable: 381-940-R



Fine Needle Electrode
 Disposable: 381-942-S Reusable: 381-942-R



Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	381-944-S	381-946-S	381-948-S	381-950-S	381-952-S
Reusable	381-944-R	381-946-R	381-948-R	381-950-R	381-952-R



Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	381-954-S	381-956-S	381-958-S	381-960-S	381-962-S
Reusable	381-954-R	381-956-R	381-958-R	381-960-R	381-962-R



Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	381-964-S	381-966-S	381-968-S	381-970-S	381-972-S	381-974-S
Reusable	381-964-R	381-966-R	381-968-R	381-970-R	381-972-R	381-974-R



Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	381-976-S	381-978-S	381-980-S	381-982-S	381-984-S	381-986-S
Reusable	381-976-R	381-978-R	381-980-R	381-982-R	381-984-R	381-986-R



Blade Electrode
Disposable: 381-1000-S Reusable: 381-1000-R



Blade Electrode Crv
Disposable: 381-1002-S Reusable: 381-1002-R



Needle Electrode
Disposable: 381-1004-S Reusable: 381-1004-R



Needle Electrode Crv
Disposable: 381-1006-S Reusable: 381-1006-R



Ball Electrode 2mm
Disposable: 381-1008-S Reusable: 381-1008-R



Ball Electrode Crv 2mm
Disposable: 381-1010-S Reusable: 381-1010-R



Ball Electrode 3mm
Disposable: 381-1012-S Reusable: 381-1012-R



Ball Electrode 3mm Crv
Disposable: 381-1014-S Reusable: 381-1014-R



Ball Electrode 4mm
Disposable: 381-1016-S Reusable: 381-1016-R



Ball Electrode 4mm Crv
Disposable: 381-1018-S Reusable: 381-1018-R



Ball Electrode 5mm
Disposable: 381-1020-S Reusable: 381-1020-R



Ball Electrode 5mm Crv
Disposable: 381-1022-S Reusable: 381-1022-R



Blade Electrode Shaft 13cm
Disposable: 381-1024-S Reusable: 381-1024-R



Needle Electrode Shaft 13cm
Disposable: 381-1026-S Reusable: 381-1026-R



Ball Electrode Shaft 3mm 13cm
Disposable: 381-1028-S Reusable: 381-1028-R



Ball Electrode 4mm Shaft 13cm
Disposable: 381-1030-S Reusable: 381-1030-R



Ball Electrode 5mm 13cm
Disposable: 381-1032-S Reusable: 381-1032-R



Circle Loop Electrode 4mm 0.7cm
Disposable: 381-1034-S Reusable: 381-1034-R



Circle Loop Electrode 6mm 0.7cm
Disposable: 381-936-S Reusable: 381-936-R



Circle Loop Electrode 8mm 0.7cm
Disposable: 381-1038-S Reusable: 381-1038-R



Circle Loop Electrode 10mm 0.7cm
Disposable: 381-940-S Reusable: 381-940-R



Plate Electrode Small
Disposable: 381-1042-S Reusable: 381-1042-R



Plate Electrode Large
Disposable: 381-1044-S Reusable: 381-1044-R



Size	1.0x0.4cm	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm
Disposable	381-1046-S	381-1048-S	381-1050-S	381-1052-S	381-1054-S
Reusable	381-1046-R	381-1048-R	381-1050-R	381-1052-R	381-1054-R



Size	1.0x0.5cm	1.0x0.7cm	1.0x0.8cm	1.0x1.0cm	1.0x1.5cm
Disposable	381-1056-S	381-1058-S	381-1060-S	381-1062-S	381-1064-S
Reusable	381-1056-R	381-1058-R	381-1060-R	381-1062-R	381-1064-R



Size	1.5x0.5cm	1.5x0.7cm	1.5x0.8cm	1.5x1.0cm	1.5x1.2cm	1.5x1.5cm
Disposable	381-1066-S	381-1068-S	381-1070-S	381-1072-S	381-1074-S	381-1076-S
Reusable	381-1066-R	381-1068-R	381-1070-R	381-1072-R	381-1074-R	381-1076-R



Size	2.0x0.8cm	2.0x1.0cm	2.0x1.2cm	2.0x1.3cm	2.0x1.5cm	2.0x2.0cm
Disposable	381-1078-S	381-1080-S	381-1082-S	381-1084-S	381-1086-S	381-1088-S
Reusable	381-1078-R	381-1080-R	381-1082-R	381-1084-R	381-1086-R	381-1088-R



Haggai Surgical Cables are available in disposable and reusable grade. We offer high quality silicon cables even for our single-use instruments thus making sure that all cables are able to withstand electrical current above the standard requirements. Furthermore, numerous tests including continuity and leak test ensures that our cables will work under the toughest requirements. Multiple sizes and plugs are available as per the need of the customer. Our cables are durable, flexible and able to withstand sterilization processes if needed. By following international standards we ensure that our cables are compatible with most major ESU's. Our cables are available in lengths 3 m, 4 m or 5 m.




Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**
 Generator Connector: **4mm Banana Plug**
 Reusable Silicon 3m Cable **HS-381-1100**
 Single Use 3m Cable **HS-381-1102**




Bipolar Cable

Instruments Connector: **European Flat Plug**
 Generator Connector: **4mm Banana Plug**
 Reusable Silicon 3m Cable **HS-381-1104**
 Single Use 3m Cable **HS-381-1106**



Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**
 Generator Connector: **Pin Distance 29mm**
 Reusable Silicon 3m Cable **HS-381-1112**
 Single Use 3m Cable **HS-381-1114**



Bipolar Cable

Instruments Connector: **European Flat Plug**
 Generator Connector: **Pin Distance 29mm**
 Reusable Silicon 3m Cable **HS-381-1108**
 Single Use 3m Cable **HS-381-1110**



Bipolar Cable

Instruments Connector: **Bipolar Artery Sealer**
 Generator Connector: **4mm Banana Plug**
 Reusable Silicon 3m Cable **HS-381-1116**



European Bipolar Cable

Instruments Connector: **European Flat Plug**
Generator Connector: **Martin Type**
Reusable Silicon 3m Cable **HS-381-1118**



US Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**
Generator Connector: **Martin Type**
Reusable Silicon 3m Cable **HS-381-1120**



European Bipolar Cable

Instruments Connector: **European Flat Plug**
Generator Connector: **Erbe Type**
Reusable Silicon 3m Cable **HS-381-1122**



US Bipolar Cable

Instruments Connector: **USA 2-Pin Plug**
Generator Connector: **Erbe Type**
Reusable Silicon 3m Cable **HS-381-1124**



Bipolar Cable

Instruments Connector: **Ethicon**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **HS-381-1126**



Monopolar Cable

Instruments Connector: **Screw Type**
Generator Connector: **Standard Bovie 8mm Plug**
Reusable Silicon 3m Cable **HS-381-1128**



Monopolar Cable

Instruments Connector: **4mm Female Fitting**
Generator Connector: **Standard Bovie 8mm Plug**
Reusable Silicon 3m Cable **HS-381-1130**



Monopolar Endoscopic Cable

Instruments Connector: **4mm Female Fitting**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **HS-381-1132**



Diathermy Cable

Instruments Connector: **4mm Male Fitting**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **HS-381-1134**



Diathermy Cable

Instruments Connector: **4.8mm Female Fitting**
Generator Connector: **4mm Banana Plug**
Reusable Silicon 3m Cable **HS-381-1136**



Diathermy Cable

Instruments Connector: **4.8mm Female Fitting**
Generator Connector: **Stan. Bovie 8mm Plug**
Reusable Silicon 3m Cable **HS-381-11HS-38**

HS-381-1150
Artery Sealer 6.5"
HS-381-1152
Artery Sealer 7.5"
HS-381-1154
Artery Sealer 8.5"
HS-381-1156
Artery Sealer 9.5"
HS-381-1158
Artery Sealer 10.5"



HS-381-1160
Non-Stick Artery Sealer 6.5"
HS-381-1162
Non-Stick Artery Sealer 7.5"
HS-381-1164
Non-Stick Artery Sealer 8.5"
HS-381-1166
Non-Stick Artery Sealer 9.5"
HS-381-1168
Non-Stick Artery Sealer 10.5"



HAGGAI



SURGICAL

IF YOU WANT TO KNOW MORE INFORMATION ABOUT OUR COMPANY, CATALOGUES, PRICES OR ANY OTHER DETAILS. PLEASE CONTACT US THROUGH THE SYSTEM.

CONTACT DETAILS



www.haggaisurgical.com



haggaisurgical@gmail.com

haggaisurgicalinternational@gmail.com



[+923246042453](https://wa.me/+923246042453)

[+923246042419](https://wa.me/+923246042419)



SADDAR CANTT DHALAY SIALKOT 51310 PUNJAB PAKISTAN